

EXHIBIT 10

BenchMarkTM XT/DiscoveryTM XT/BenchMark LT Service Manual

**CONFIDENTIAL
ATTORNEYS ONLY**

VM 000425

BenchMark™ XT/Discovery™ XT/ BenchMark LT Service Manual

This manual's Part Number is 2000300 Revision B.

**CONFIDENTIAL
ATTORNEYS ONLY**

VM 000426

Copyright © Ventana Medical Systems, Inc. 2004
All Rights Reserved.

**CONFIDENTIAL
ATTORNEYS ONLY**

VM 000427

BenchMark XT, LT/Discovery XT Service Manual

2.0 WHAT BENCHMARK XT/ DISCOVERY XT/ BENCHMARK LT DOES

2.1 Tissue Staining

The BenchMark XT/Discovery XT/BenchMark LT instrument automates an exacting and somewhat tedious laboratory procedure used to produce visible *stains* on biological tissue samples.

- These stained tissue samples (on glass slides) are used by a pathologist to diagnose various diseases
- The chemistries involved are called immunohistochemistry (IHC) or *in situ* hybridization (ISH)

2.2 Tissue Staining Process

A multi-step process is needed to produce a stain. In each step, a chemical called a *reagent* is applied to the tissue sample.

- The reagent and the tissue sample undergo a biochemical reaction
 - ♦ The biochemical reaction must be allowed to occur over a specified amount of time
 - ♦ The biochemical reaction must be allowed to occur at a precisely controlled temperature
- The process is then repeated with another reagent, which requires yet another reaction time
- A number of steps spaced over one to four hours is usually required to produce the desired stain

The BenchMark XT/Discovery XT instrument allows up to 30 slides to be loaded and stained in a single run, while the BenchMark LT instrument allows up to 20 slides to be loaded and stained in a single run. All three instruments can use up to 35 different reagents. A laboratory technician can specify the stain needed for each slide, load the reagents, then start the instrument and simply walk away until all of the stains have been completed.

BenchMark XT, LT/Discovery XT Service Manual

- This is called **bar code blow-off**
 - ◇ Liquid on the bar code label prevents the bar code reader from reading
- ◆ Create an LCS puddle separated from the bar code label
 - A jet of air called an **air knife** keeps the puddle from leaking or wicking into the bar code label
 - ◇ Once separated, the edge of the bar code label will act as a dam, forcing the puddle to remain intact on the slide like a drop of water on a window
- **DISPENSING**—The specified reagent (identified by its bar code label) is dispensed into the puddle created in the previous steps
 - ◆ The LCS solution floats to the top, again covering the solution, which is now a mixture of buffer and reagent
- **MIXING**—The instrument uses **vortex mixers** (air streams) to mix the solution and distribute it over the tissue sample on the slide
 - ◆ The mixing action is performed several times
 - ◆ The mixing action helps the biochemical reaction
 - ◆ The mixing action ensures that the reagents come into contact with the entire tissue sample on the slide

Temperature is an important ingredient in the process.

- BenchMark XT/Discovery XT/BenchMark LT has slide ThermoPads™ to control the temperature of individual slides so that the chemical reactions will occur under optimum conditions

During a *run* of the instrument, the six step sequence above is repeated according to the protocol specified for *each* slide.

- Reagent selection and the time between dispense steps are part of the protocol for a given slide
- Each slide the BenchMark XT/Discovery XT/ BenchMark LT instrument can hold may have a separate protocol, or all may have the same protocol
- BenchMark XT/Discovery XT/BenchMark LT obviously saves a lot of laboratory manpower and potential for error compared to doing it all by hand

BenchMark XT, LT/Discovery XT Service Manual

3.0 MAJOR SYSTEM MODULES

BenchMark XT/Discovery XT/BenchMark LT is comprised of four major modules:

- The stainer subassembly
- The automated fluidics subassembly
- The waste subassembly
- The personal computer

3.1 Stainer Subassembly

This subassembly is where the chemistry is performed.

- The BenchMark XT and Discovery XT can hold up to 30 slides
- The BenchMark LT can hold up to 20 slides
 - ♦ Each of the slide positions has a fixed integrated slide holder and ThermoPad
- It has a rotating nozzle plate that is rotated by a stepper motor
 - ♦ The rotary position of the nozzle plate is relative to a home position obtained from a sensor
 - ♦ The nozzle plate incorporates
 - A slide bar code reader
 - A variety of valves and nozzles for dispensing and removing fluids
- It has a reagent carousel with spaces for 35 reagent dispensers
 - ♦ The reagent carousel is rotated by a stepper motor
 - The rotary position of the reagent carousel is relative to a home position obtained from a sensor
 - ♦ Reagents are dispensed from disposable containers that snap into the **dispenser plate** on the lower part of the reagent carousel
 - To dispense a reagent, BenchMark XT/Discovery XT/BenchMark LT rotates the reagent into position below a plunger which presses down on a syringe built into the container

Staining Module

- Reagent Dispensers
- Bar Code Readers
- Liquid and Air Nozzles
 - ♦ Nozzle Arrangement
 - ♦ Dual Rinse Nozzles
 - ♦ Jet Drain Nozzle
 - ♦ Multispense Nozzles
 - ♦ LCS Nozzles
 - ♦ Air Knife Nozzle
 - ♦ Vortex Mixer Nozzles
- Restrictors and Filters
- Fan and Sensor
- Slide Tray Sensor
- Tub and Overflow Sensor
- Instrument Control Panel
- Power Supply
- Cables and Wiring Harnesses
- Plumbing

4.1 Reagent Carousel and Nozzle Plate

BenchMark XT/Discovery XT/BenchMark LT uses a rotating carousel to carry reagent dispensers. It also uses a rotating nozzle plate for dispensing and removing fluids, and reading slide bar codes.

4.1.1 Reagent Carousel

The carousel is able to rotate either clockwise or counterclockwise.

- The reagent carousel is rotated by a power transmission assembly
 - ♦ It is positioned by a stepper motor through a belt drive and sprockets

A dispenser plate holds the reagent dispensers.

- The microcontroller is able to operate the stepper in either a low-torque or a high-torque mode

Staining Module




Figure 14. Outer Mixer

- The inner and outer vortex mixers create a gentle vortex effect over the liquid puddle on the slide
 - ♦ As the nozzle plate rotates, the vortex mixers swirl the liquid puddle on the slide clockwise, then counterclockwise at several stations
 - ♦ This stirring action
 - Thoroughly mixes the reagents and antibodies
 - Provides uniform wetting of the specimen during incubation